



### Advance Electric Table (UMS - E-201)

- ✓ SS-304
- ✓ Surgeons Panel
- ✓ Electric (remote operated)
- ✓ C Arm Compatible
- ✓ Floor Locking
- ✓ Zero Positioning

#### Features :-

- ☞ Fully electrically operated table, five sections, radio-transparent table top with eccentric pillars and Die Compressed base.
- ☞ Interchangeable head and leg sections (manually).
- ☞ All positions operated through cable remote.
- ☞ Electric floor locking.
- ☞ Inbuilt kidney bridge
- ☞ Auto leveling (zero function)
- ☞ We are providing two Handsets. Additional fixed control panel provided on the table top, which can be used in case of failure of main hand set.
- ☞ Indication of Battery Back provided in th hand set.
- ☞ Handset is providing Dual Memory Function which can be set and reset.

#### TECHNICAL DATA

Top dimension	L 1900mm x W 522mm
Height adjustment	750 mm - 1000 mm
Trendelenburg / Reverse	30° / 25°
Lateral till	20° / 20°
Kidney elevator	130 mm
Back Rest ( up/down)	60° / 25°
Leg Rest ( up/down)	10° / 90°
Head Rest ( up/down)	20° / 60°
Power Supply	24V DC
Battery Backup	2-3 hours (optional)
Mains voltage	AC 220V 10% 50Hz
Patient Weight Capacity	250 kgs