



Advance Electric Table (UMS - E-201)

- ✓ SS-304
- ✓ Surgeons Panel
- ✓ Electric (remote operated)
- ✓ C Arm Compatible
- ✓ Floor Locking
- ✓ Zero Positioning

Features :-

- ☞ Fully electrically operated table, five sections, radio-transparent table top with eccentric pillars and Die Compressed base.
- ☞ Interchangeable head and leg sections (manually).
- ☞ All positions operated through cable remote.
- ☞ Electric floor locking.
- ☞ Inbuilt kidney bridge
- ☞ Auto leveling (zero function)
- ☞ We are providing two Handsets. Additional fixed control panel provided on the table top, which can be used in case of failure of main hand set.
- ☞ Indication of Battery Back provided in th hand set.
- ☞ Handset is providing Dual Memory Function which can be set and reset.

TECHNICAL DATA

Top dimension	L 1900mm x W 522mm
Height adjustment	750 mm - 1000 mm
Trendelenburg / Reverse	30° / 25°
Lateral till	20° / 20°
Kidney elevator	130 mm
Back Rest (up/down)	60° / 25°
Leg Rest (up/down)	10° / 90°
Head Rest (up/down)	20° / 60°
Power Supply	24V DC
Battery Backup	2-3 hours (optional)
Mains voltage	AC 220V 10% 50Hz
Patient Weight Capacity	250 kgs



Hydraulic C-ARM Table (UMS - CH02)

- ✓ SS-304
- ✓ Surgeons Panel
- ✓ Hydraulic Operated
- ✓ C Arm Compatible

Features :-

- ☞ Hydraulic and Mechanical operated table, five sections, radio-translucent table top pillars and Die Compressed base.
- ☞ The base is eccentric for trouble-free C-ARM movements.
- ☞ Leg and Head section are interchangeable.
- ☞ Eccentric pillars, Stainless Steel sheet fitting.
- ☞ Detachable head, leg and pelvic section.
- ☞ Inbuilt kidney bridge.
- ☞ Five sectional radio-translucent.
- ☞ Base and cover made of attractive, easy to clean stainless steel.
- ☞ Sophisticated mechanics provide smooth, stepless articulation of the table top for precise patient positioning.

TECHNICAL DATA

Top dimension	L 1900mm x W 522mm
Height adjustment	750 mm - 1000 mm
Trendelenburg / Reverse	30° / 25°
Lateral tilt	20° / 20°
Kidney elevator	130 mm
Back Rest (up/down)	60° / 25°
Leg Rest (up/down)	10° / 90°
Head Rest (up/down)	20° / 60°
Power Supply	24V DC
Battery Backup	2-3 hours (optional)
Mains voltage	AC 220V 10% 50Hz
Patient Weight Capacity	250 kgs



Electronic Surgical Operation Table (UMS - ES01)

Electrically operated table, five sections, stainless steel table top having all surgical position i.e. treandelenburg, reverse treandelenburg, lateral tillt, flex & chair position, interchangeable head and leg sections. All positions operated through cranks fitted at the side end of the table. And raising / lowering through cable remote. Supplied complete with standard ss accessories & matterss.

TECHNICAL DATA

Top dimension	L 1900mm x W 522mm
Height adjustment	750 mm - 1000 mm
Trendelenburg / Reverse	30° / 25°
Lateral till	20° / 20°
Kidney elevator	130 mm
Back Rest (up/down)	60° / 25°
Leg Rest (up/down)	10° / 90°
Head Rest (up/down)	20° / 60°
Power Supply	24V DC
Battery Backup	2-3 hours (optional)
Mains voltage	AC 220V 10% 50Hz
Patient Weight Capacity	250 kgs



Hydraulic Surgical Operation Table (UMS - HS02)

Hydraulic operated table, five sections, stainless steel table top having all surgical position i.e. raising/lowering, treandelenburg, reverse treandelenburg, lateral tillt, flex/re-flex & chair position, interchangeable head and leg sections. All positions operated through cranks fitted under the head end of the table. In build kidney bridge Supplied complete with standard ss accessories & matterss.



**Unicare Remote Oprated
LED Light Round (UMS - RF1)**

Luminous power 175000 lux Single Dome
Color temp 5500k,
Illumines field 15 cm
Color rendering index RS-95
Electronic transformer
Supply voltage : 230v 50Hz AC
Lamp housing polycarbontate
Class of protection : 1CE
Cold light & white light
No black spot in the center
Superior lamp life of 55000hr,
or 2 year manufacturing warranty.



**Unicare Remote Oprated
LED Light Apple (UMS - RF1)**

Luminous power 16000 lux Single
Dome
Color temp 5500k,
Illumines field 15 cm
Color rendering index RS-95
Electronic transformer
Supply voltage : 230v 50Hz AC
Lamp housing polycarbontate
Class of protection : 1CE
Cold light & white light
No black spot in the center
Superior lamp life of 55000hr,
or 2 year manufacturing warranty.



**HEX LED OT Light Double Dome
(UMS - HL8)**

Luminous power - 160000 + lux Single
Dome
Color temp 5500K
Illumines field 15
Color rendering index RS-95
Electronic transformer
Supply voltage : 230V 50hz AC
Lamp housing polycarbonate
Class of protection : 1 CE
Cold light & white light
No black spot in the center
Superior lamp life of 35000hr
or 2 year manufacturing warranty.