



Attandat Bed cum chair

**Technical Specifications**

Overall Approximate Dimension	
Dimension From End to End*	1850 mm x 780 mm, 73" x 31" (L x W)
Platform Dimension*(Without Mattress)	1850 mm x 640 mm, 73" x 25" (L x W)
Fixed Height*	515 mm, 20" H (without upholstery)
Upholstery Thickness*	75 mm, 3" H
Backrest Angle**	0° - 70°
Safe Working Load	135 Kgs
Patient Load Bearing Capacity	135 Kgs