



**ATTANDANT BED**

Overall size: 1800L X 600W X 450H mm CRCA ms tube frame work  
 With perforated CRCA sheet finis pre-treated epoxy powder coated.



**FOLDING ATTENDANT BED**

Overall size: 1800L X 600W X 450H mm two part folding adjustment ratchet through CRCA ms tube frame work With perforated CRCA sheet Finis pre-treated epoxy powder coated.



**Attendant Stool (UMS-536)**



**Revolving Stool (UMS-P537)**



**Revolving Stool SS (UMS-537)**



**Bed Side Screen (UMS-404)**



**Slime Stand (UMS-S405)**