



**ELECTRONIC CENTRAL LOCKIN (5) FUCTION ICU BED UMS-ED4**

Technical features		
Bed structure material	MS tube 25x75 pipe 1.4 mm thickness and 1.25 mm thickness CRCA sheet with perforated copper wire c02 welding	(1)Head position up down(2) foot position up down(3) height hi lo position (4) UP tredelenburg angle(5) reverse tredelenburg angle
Structure material	MS tube pipe 1.4 mm and 1.25m CRCA sheet	ABS Moulded Large size panels
Overall bed outer damnation	2100L X980W	ABS 4pis large size tuck away safety side railing
Mattress platform damnation	2000L X 880W	125 mm catore wheel with central locking systems
Height adjustment	450 to 750 mm electronic linear actuators through	CRCA perforated sheet
Reveres trengle	Electronic linear actuator through	Electronic actuator
Up trengle	Electronic linear actuator through	Central locking systems
Operating systems	Electronic actuator By electronic hand remote	4 side bed safety side wafer
Product finis	Pre-treated epoxy powder coated	4 side iv pole hole
Safe working load	200 kg	Mattress IV stand cylinder trolley optional