



**MONITOR TROLLEY UMS-512**  
 Overall size 700L X 500D X1580H  
 Adjustable 4 shelf between  
 Top shelf and drawer drawer are  
 available with lock and key  
 mechanism  
 Four 100mm non rusted castor  
 wheel with 2 brak  
 Finis pre-treated and epoxy  
 powder coated also available  
 complete steel



**Crash Cart Trolley Ums-513**  
 Overall Approx .size=1000mm L X 510mm D  
 X 1540mm H  
 SS Tubular Frame Work  
 Ss Top, Middle And Bottom Tray  
 Six Removal Bins  
 Two Alkon Make Modular Storage Units With  
 Lock Three Drawers Each  
 Oxygen Cylinder Holder  
 Ss Iv Rod And Cardiacmassage Board  
 Trolley Mountedon Four Non- Rusting 100  
 Mm Dia Non Rusting Wheel 2 With Brakes  
 And 2 Without Brakes  
 Alternative 513p In Expoy Powder Coated  
 Framework And Trays



**EMERGENCY MEDICINE  
 TROLLEY COMPLETE SS  
 UMS- S514**  
 Overall Approximate  
 SIZE.=990MM L X 490MM D  
 X 1600MM H  
 SS TUBULAR  
 FRAMEWORK  
 6 LARG SLIDING DRAWER  
 one side fitted plastic bucket  
 SS IV Road easy mobavel  
 handle and trolley mounted  
 on on four non rusting  
 100mm castor wheel  
 With 2 brack



**EMERGENCY MEDICINE  
 TROLLEY POWDER COATED  
 UMS-P514**  
 Overall Approximate SIZE.=990MM  
 L X 490MM D X 1600MM H  
 CRCA tubular frame work 3 drawer  
 medium and 2 larg size slidind  
 drawerSS IV Road easy movable  
 SS handle and trolley mounted on  
 on four non rusting 100mm castor  
 wheel  
 With 2 brack