



DELUX EXAMINATION COUCH UMS-D500

Technical specification & Features

*Overall approximate dimension 1940L x 600w x 850H mm CRCA rectangular tube frame with CRCA sheet cabinets Fixed upholstered top in two section Retractable with 13 unique design drawers Headrest adjustable on gas spring Three cabinets and three drawers B.P apparatus tray pre treated and powder coated



EXAMINATION COUCH UMS-502

Features

Overall approximate dimension 1940L x 600w x 850H mm CRCA rectangular tube frame with CRCA sheet cabinets Fixed upholstered top in two section With U cut Retractable step stool Head rest adjustable on gas spring Three cabinets and three drawers BP apparatus tray pre treated and powder coated



GYNIC EXAMINATION COUCH UMS-501

Technical specification Features

*Overall approximate dimension 1940L x 600w x 850H mm CRCA rectangular tube frame with CRCA sheet cabinets Fixed upholstered top in two section Retractable step stool Headrest adjustable on gas spring Three cabinets and three drawers BP apparatus tray pre treated and powder coated