



Also Available MS top



ABS with Transparent Top



X-Ray Cassette



Hi-Lo Position



Side Railing up down

ATTENDANT CONTROL PANEL

ES-12DW




Weighing & Monitoring Module integrated in a complete 5 function ACP.

ES-12HD



Complete 5 function ACP which including Lock-Out, Auto-Regression controls.



ALIXIGER ELECTRONIC TILT POSITION ICU BED (7) FUCTION UMS-PT01

Technical features		STANDARD CONFIGURATION
Bed structure material	MS tube 25x75 pipe 1.4 mm thickness and 1.25 mm thickness CRC sheet with perforated copper wire c02 welding	(1)Head position up down(2) foot position up down(3) height HI-LO position (4) UP tredelenburg angle(5) reverse tredelenburg angle(6)lateral tillt position
Structure material	MS tube pipe 1.4 mm and 1.25 mm CRC sheet	ABS Moulded large panels With nursing operating system
waiss structure	ABS moulded cover fitted	
Overall bed outer damnation	2200L X1020W	ABS 4pis large size tuck away safety side railing With patient operating remote
Mattress platform damnation	2000L X 900W	125 mm catore wheel with central locking systems
Height adjustment	400 to 750 mm by linear actuators	ABS moulded platform and lower bed fram ABS cover
Reveres trengleburg	By linear actuators	Electronic actuator
Up trengleburg	By linear actuators	Electronic Central locking systems
Operating systems	By electronic cylindrical hand remote	4 side bed safety side wafer
Product finis	Pre-treated and epoxy powder coated	4 side iv pole hole
Safe working load	300 kg	Mattress IV stand cylinder trolley