



### Hydraulic C-ARM Table (UMS - CH02)

- ✓ SS-304
- ✓ Surgeons Panel
- ✓ Hydraulic Operated
- ✓ C Arm Compatible

#### Features :-

- ☞ Hydraulic and Mechanical operated table, five sections, radio-translucent table top pillars and Die Compressed base.
- ☞ The base is eccentric for trouble-free C-ARM movements.
- ☞ Leg and Head section are interchangeable.
- ☞ Eccentric pillars, Stainless Steel sheet fitting.
- ☞ Detachable head, leg and pelvic section.
- ☞ Inbuilt kidney bridge.
- ☞ Five sectional radio-translucent.
- ☞ Base and cover made of attractive, easy to clean stainless steel.
- ☞ Sophisticated mechanics provide smooth, stepless articulation of the table top for precise patient positioning.

#### TECHNICAL DATA

Top dimension	L 1900mm x W 522mm
Height adjustment	750 mm - 1000 mm
Trendelenburg / Reverse	30° / 25°
Lateral tilt	20° / 20°
Kidney elevator	130 mm
Back Rest ( up/down)	60° / 25°
Leg Rest ( up/down)	10° / 90°
Head Rest ( up/down)	20° / 60°
Power Supply	24V DC
Battery Backup	2-3 hours (optional)
Mains voltage	AC 220V 10% 50Hz
Patient Weight Capacity	250 kgs