



**MONITOR TROLLEY UMS-512**  
 Overall size 700L X 500D X1580H  
 Adjustable 4 shelf between  
 Top shelf and drawer drawer are  
 available with lock and key  
 mechanism  
 Four 100mm non rusted castor  
 wheel with 2 brak  
 Finis pre-treated and epoxy  
 powder coated also available  
 complete steel



**Crash Cart Trolley Ums-513**  
 Overall Approx .size=1000mm L X 510mm D  
 X 1540mm H  
 SS Tubular Frame Work  
 Ss Top, Middle And Bottom Tray  
 Six Removal Bins  
 Two Alkon Make Modular Storage Units With  
 Lock Three Drawers Each  
 Oxygen Cylinder Holder  
 Ss Iv Rod And Cardiacmassage Board  
 Trolley Mountedon Four Non- Rusting 100  
 Mm Dia Non Rusting Wheel 2 With Brakes  
 And 2 Without Brakes  
 Alternative 513p In Expoy Powder Coated  
 Framework And Trays



**EMERGENCY MEDICINE  
 TROLLEY COMPLETE SS  
 UMS- S514**  
 Overall Approximate  
 SIZE.=990MM L X 490MM D  
 X 1600MM H  
 SS TUBULAR  
 FRAMEWORK  
 6 LARG SLIDING DRAWER  
 one side fitted plastic bucket  
 SS IV Road easy mobavel  
 SS handle and trolley mounted  
 on on four non rusting  
 100mm castor wheel  
 With 2 brak



**EMERGENCY MEDICINE  
 TROLLEY POWDER COATED  
 UMS-P514**  
 Overall Approximate SIZE.=990MM  
 L X 490MM D X 1600MM H  
 CRCA tubular frame work 3 drawer  
 medium and 2 larg size slidind  
 drawerSS IV Road easy movable  
 SS handle and trolley mounted on  
 on four non rusting 100mm castor  
 wheel  
 With 2 brak



**INSTRUMENT CABINET WITH DRAWER BOX UMS-515**

Overall size 1600H X 900W X 500D mm top side 2 self with glass door  
Midil position 2 sliding draz lower machine pressed box with 2 door easy  
movable SS handle complete stainless steel also available pre treated  
powder coated finis



**INSTRUMENT CABINET DUBLE DOOR UMS-516**

Overall size 1600H X 900W X 500D mm top side 3 adjustable self  
with glass fitted duble door available with hand lock key mechanism  
SS handle complete stainless steel also available pre treated powder coated  
finis



**INSTRUMENT TROLLEY (2) SELVES UMS-517**

(A) overall size : 27L X 18W X 36H inch  
 (B) overall size: 30L X 24W X 36H inch  
 (C) overall size: 42L X 24W X 36H inch  
 stainless steel tubular frame mounted on 100 mm non rusted castor wheel with 2 moulded SS self with protective 3 side railings fitted



**INSTRUMENT TROLLEY (3) SELVES UMS-518**

(A) overall size: 30L X 24W X 36H inch  
 (B) overall size: 42L X 24W X 42H inch  
 stainless steel tubular frame mounted on 100 mm non rusted castor wheel with 3 moulded SS self with protective 3 side railings fitted



**DRESSING TROLLEY (2) SELVES UMS-519**

(A) overall size : 27L X 18W X 36H inch  
 (B) overall size: 30L X 24W X 36H inch  
 stainless steel tubular frame mounted on 100 mm non rusted castor wheel with 2 moulded SS self with protective 3 side railing bowl bucket fitted



**OT LINEN TROLLEY UMS-520**

(A) overall size : 36L X 18W X 36H inch easy movable handle stainless steel tubular frame mounted on 100 mm non rusted castor wheel with protective 4 side pipe frame fitted canvas cover



**MAYO INSTRUMENT TROLLEY DUBLE PILLAR UMS-521**

Overall size: tray size 36L X 20W 30 TO 50 inch adjustable height  
 SS tubular frame mounted on four non rusted castor wheel with 2 brake



**MAYO INSTRUMENT TROLLEY SINGLE PILLAR UM-522**

Overall size: tray size 22L X 16W 30 TO 50 inch adjustable height  
 SS tubular frame mounted on four non rusted castor wheel with 2 brake