



**Unicare Remote Oprated  
LED Light Round (UMS - RF1)**

Luminous power 175000 lux Single Dome  
Color temp 5500k,  
Illumines field 15 cm  
Color rendering index RS-95  
Electronic transformer  
Supply voltage : 230v 50Hz AC  
Lamp housing polycarbontate  
Class of protection : 1CE  
Cold light & white light  
No black spot in the center  
Superior lamp life of 55000hr,  
or 2 year manufacturing warranty.



**Unicare Remote Oprated  
LED Light Apple (UMS - RF1)**

Luminous power 16000 lux Single  
Dome  
Color temp 5500k,  
Illumines field 15 cm  
Color rendering index RS-95  
Electronic transformer  
Supply voltage : 230v 50Hz AC  
Lamp housing polycarbontate  
Class of protection : 1CE  
Cold light & white light  
No black spot in the center  
Superior lamp life of 55000hr,  
or 2 year manufacturing warranty.



**HEX LED OT Light Double Dome  
(UMS - HL8)**

Luminous power - 160000 + lux Single  
Dome  
Color temp 5500K  
Illumines field 15  
Color rendering index RS-95  
Electronic transformer  
Supply voltage : 230V 50hz AC  
Lamp housing polycarbonate  
Class of protection : 1 CE  
Cold light & white light  
No black spot in the center  
Superior lamp life of 35000hr  
or 2 year manufacturing warranty.