



PEDIATRIC BED DELUX UMS-P01

- Overall approximate dimension: 1490 mm x 830 mm x 600 mm, 59" x 33" x 24" (L x W x H).
- Bed frame approximate dimension: 1370 mm x 760 mm, 54" x 30" (L x W).
- Perforated M.S. sheet top.
- Both bows of 1067 mm (42") H, with one horizontal, and seven vertical supports and rubber feet.
- Full length drop side rails.
- Provision for mosquito net poles.
- Pre-treated and powder coated.
- Supplied in KDC.

Optional Accessories

Rexine covered mattress of 100 mm thickness (0111D)



PEDIATRIC BED DELUX UMS-P02



BABY CRADEL UMS B03

- Overall approximate dimension: 920 mm x 455mm x 950 mm, 36" x 18" x 37" (L x W x H).
- Approximate crib size: 760 mm x 380 mm x 305 mm, 30" x 15" x 12" (L x W x H).
- 3" PU Foam Mattress with PVC Rexine.
- Tubular frame fitted with swiveling castors of 125 mm diameter synthetic body castors two with brake and two without brake.
- Pre-treated and epoxy powder coated.