



FUTURA ELECTRIC STD (2) FUNCTION BEDUMS-ED7

| Technical features | | Standard configuration |
|-----------------------------|---|---|
| Bed structure material | MS tube 25x75 pipe 1.4 mm thickness and 1.25 mm thickness CRC sheet with perforated copper wire c02 welding | (1)Head position up down (2) foot position up down |
| Structure material | MS tube pipe 1.4 mm and 1.25 mm CRC sheet | ABS Moulded panels |
| Overall bed outer damnation | 2000L X980W | Aluminium collapsible safety side railing |
| Mattress platform damnation | 1900L X 880W | 100 mm castor wheel with 2 brake |
| Height adjustment | Fix height 600 mm | MS CRC perforated sheet |
| | | Electronic actuator |
| | | 4 wheel With 2 brake systems |
| Operating systems | Electronic actuator By electronic hand remote | 4 side bed safety side wafer |
| Product finis | Pre-treated and epoxy powder coated | 4 side iv pole hole |
| Safe working load | 200 kg | Mattress IV stand cylinder trolley optional |