



FUTURA ELECTRIC SEMI ICU (2) FUCTION ICU BED-ED6

Technical features		Standard configuration
Bed structure material	MS tube 25x75 pipe 1.4 mm thickness and 1.25 mm thickness CRC sheet with perforated copper wire c02 welding	(1)Head position up down by linear electronic actuators (2) foot position up down by linear electronic actuators
Structure material	MS tube pipe 1.4 mm and 1.25 mm CRC sheet	ABS Moulded panels
Overall bed outer damnation	2100L X980W	ABS 4pis large size tuck away safety side railing
Mattress platform damnation	2000L X 880W	125 mm castor wheel with 2 brake
Height adjustment	Fix height 600 mm	MS CRC perforated sheet
		Electronic actuator
		4 wheel With 2 brake systems
Operating systems	Electronic actuator By electronic hand remote	4 side bed safety side wafer
Product finis	Pre-treated epoxy powder coated	4 side iv pole hole
Safe working load	200 kg	Mattress IV stand cylinder trolley optional