



DELUX EXAMINATION COUCH UMS-D500

Technical specification & Features

*Overall approximate dimension 1940L x 600w x 850H mm CRCA rectangular tube frame with CRCA sheet cabinets Fixed upholstered top in two section Retractable with 13 unique design drawers Headrest adjustable on gas spring Three cabinets and three drawers B.P apparatus tray pre treated and powder coated



EXAMINATION COUCH UMS-502

Features

Overall approximate dimension 1940L x 600w x 850H mm CRCA rectangular tube frame with CRCA sheet cabinets Fixed upholstered top in two section With U cut Retractable step stool Head rest adjustable on gas spring Three cabinets and three drawers BP apparatus tray pre treated and powder coated



GYNIC EXAMINATION COUCH UMS-501

Technical specification Features

*Overall approximate dimension 1940L x 600w x 850H mm CRCA rectangular tube frame with CRCA sheet cabinets Fixed upholstered top in two section Retractable step stool Headrest adjustable on gas spring Three cabinets and three drawers BP apparatus tray pre treated and powder coated



ELECTRONIC EXAMINATION TABLE UMS-EX1
 Overall approximate dimension 1940L x 600w x 600 - 800HI-LO
 Position by linear electronic actuators throw adjustable easy accessibility and paper roll holder at foot and of the couch head rest on gas spring hand rest side support 100 mm non rested castor wheel with 2 brack finis pre-treated and epoxy powder coated



EXAMINATION TABLE UMS-ET504
 Overall approximate dimension 1800L X 600W X800H Fix height
 head rest on manually ratchet throw finis pre-treated and epoxy powder coated



U CUT EXAMINATION TABLE UMS-GT503
 Overall approximate dimension 1800L X 600W X800H Fix height with U CUT head rest on manually ratchet throw finis pre-treated and epoxy powder coated



Foot Stap Delux ABS top (UMS - F2)



Foot Stap (UMS - F2)



Foot Stap Single (UMS - F1)