



ELECTRONIC DY LISYS CHAIR UMS-DC1

- Chair frame with Epoxy powder coating
- Anti-bacterial MS Frame
- Quick Positioning possible due to hand control and 3 motors
- Raising and Lowering Legrest movement
- Detachable Leg support
- Chair on 4 wheels
- Center Locking with direction locking system

Technical Data :

Length 1710 mm
Width 610 mm

Backrest Angle(Up/Down)	0°-80°
Power Supply	230 V AC 50Hz
Operating voltage of Actuators	24V DC
Operating system By Remote	

Dialysis chair

- Chair frame with Epoxy powder coating
- Anti-bacterial MS Frame
- Quick Positioning possible due to hand control and 3 motors
- Raising and Lowering Legrest movement
- Detachable Leg support
- Chair on 4 wheels
- Center Locking with direction locking system

Technical Data :

Length	1710 mm
Width	610 mm
Backrest Angle(Up/Down)	0°-80°
Manually operated Gas shocker Throw	