



SENTRA DELUX MANUALLY SEMI FOWLER BED ABS PANEL UMS-107

Technical features		Standard configuration
Bed structure material	MS tube 25x75 pipe 1.4 mm thickness and 1.25 mm thickness CRCA sheet with perforated copper wire c02 welding	(1)Head position up down
Structure material	MS tube pipe 1.4 mm and 1.25 mm CRC sheet	ABS Moulded panels
Overall bed outer damnation	2000L X980W	Aluminium collapsible safety side railing optional
Mattress platform damnation	1900L X 880W	
Height adjustment	Fix height 600 mm	MS CRC perforated sheet
		Hand key with polymer handle
Operating systems	Hand key	4 side bed safety side wafer
Product finis	Pre-treated and epoxy powder coated	4 side iv pole hole
Safe working load	200 kg	Mattress IV stand cylinder trolley

SENTRA MANUALLY SEMI FOWLER BED SS PANEL UMS-108

Technical features		Standard configuration
Bed structure material	MS tube 25x75 pipe 1.4 mm thickness and 1.25 mm thickness CRC sheet with perforated copper wire c02 welding	(1)Head position up down
Structure material	MS tube pipe 1.4 mm and 1.25 mm CRC sheet	Steel deluxe panels
Overall bed outer damnation	2000L X980W	Aluminium collapsible safety side railing optional
Mattress platform damnation	1900L X 880W	
Height adjustment	Fix height 600 mm	MS CRC perforated sheet
		Hand key with polymer handle
Operating systems	Hand key	4 side bed safety side wafer
Product finis	Pre-treated and epoxy powder coated	4 side iv pole hole
Safe working load	200 kg	Mattress IV stand t optional