



Technical Specifications

Overall Approximate Dimension	
Dimension From End to End*	1850 mm x 780 mm, 73" x 31" (L x W)
Platform Dimension*(Without Mattress)	1850 mm x 640 mm, 73" x 25" (L x W)
Fixed Height*	515 mm, 20" H (without upholstery)
Upholstery Thickness*	75 mm, 3" H
Backrest Angle**	0° - 70°
Safe Working Load	135 Kgs
Patient Load Bearing Capacity	135 Kgs



ATTANDANT BED

Overall size: 1800L X 600W X 450H mm CRCA ms tube frame work
With perforated CRCA sheet finis pre-treated epoxy powder coated.



FOLDING ATTENDANT BED

Overall size: 1800L X 600W X 450H mm two part folding adjustment ratchet through CRCA ms tube frame work With perforated CRCA sheet Finis pre-treated epoxy powder coated.



**Attendant Stool
(UMS-536)**



**Revolving Stool
(UMS-P537)**



**Revolving Stool SS
(UMS-537)**



**Bed Side Screen
(UMS-404)**



**Sline Stand
(UMS-S405)**