



Electronic Surgical Operation Table (UMS - ES01)

Electrically operated table, five sections, stainless steel table top having all surgical position i.e. trendelenburg, reverse trendelenburg, lateral tilt, flex & chair position, interchangeable head and leg sections. All positions operated through cranks fitted at the side end of the table. And raising / lowering through cable remote. Supplied complete with standard ss accessories & matterss.

TECHNICAL DATA

Top dimension	L 1900mm x W 522mm
Height adjustment	750 mm - 1000 mm
Trendelenburg / Reverse	30° / 25°
Lateral tilt	20° / 20°
Kidney elevator	130 mm
Back Rest (up/down)	60° / 25°
Leg Rest (up/down)	10° / 90°
Head Rest (up/down)	20° / 60°
Power Supply	24V DC
Battery Backup	2-3 hours (optional)
Mains voltage	AC 220V 10% 50Hz
Patient Weight Capacity	250 kgs



Hydraulic Surgical Operation Table (UMS - HS02)

Hydraulic operated table, five sections, stainless steel table top having all surgical position i.e. raising/lowering, trendelenburg, reverse trendelenburg, lateral tilt, flex/re-flex & chair position, interchangeable head and leg sections. All positions operated through cranks fitted under the head end of the table. In build kidney bridge Supplied complete with standard ss accessories & matterss.