



HYDRAULIC HEIGHT ADJUSTABLE LABUR TABLE UMS-HD03

Overall platform 1890MM X 900 X 800mm mattress platform 1800L X 650W three section MS tubular frame with SS Top Backrest adjustable on ratchet=750 mm x 627mm Hydraulic operated height adjustment 600 to 800mm adjustable height Middle section with "u" cut= 350mm x 655mm Leg section can be telescoped under the main frame completely for lithotomic positions Hand operated gas spring assisted trendilenburg position. Complete with stainless stell head and side rails hand grips, a pair of rubber padded stainless steel lithotomy rods with rexine ankle straps and stainless steel bowl under 'u' cut. Size of railing=195mm x 685 mm. Supplied in skd condition. Optional accessories Lithotomy cruches in lieu of lithotomy ro Rexine covered mattress of 25 mm Stainless steel I.V.rod – hook



OBSTETRIC LABOUR TABLE FULLY S.S. UMS-OLT04

Overall platform 1890MM X 900 X 800mm mattress platform 1800L X 650W three section SS tubular frame Backrest adjustable on ratchet=750 mm x 627mm Middle section with "u" cut= 350mm x 655mm Leg section can be telescoped under the main frame completely for lithotomic positions Hand operated gas spring assisted trendilenburg position. Complete with stainless stell head and side rails hand grips , a pair of rubber padded stainless steel lithotomy rods with rexine ankle straps and stainless steel bowl under 'u' cut. Size of railing=195mm x 685 mm. Supplied in skd condition. Optional accessoriesLithotomy cruches in lieu of lithotomy ro Rexine covered mattress of 25 mm Stainless steel I.V.rod – hook