



ABS MOULDED BED SIDE LOCKER UMS-AB02

Overall box size : 425L x 425W X 800H mm and drawer size 125HX 300W X 350D mm

Locker made from ABS moulded all side material having one box and one drawer fitted with ABS moulded top with raised edges, drawer and cabinet door front side buffer at rear to avoid collision Base fitted with 4 swivel castors of 50mm diameter made from synthetic material without brake



PRIME BED SIDE LOCKER WITH OVERBED TABLE UMS-BO01

Overall box size : 450L x 450W X 800H mm and drawer size 125HX 300W X 350D mm

Locker made from machine pressed CRCA ms sheet enclosed on three side which is

Pre-treated and epoxy powder coated finish with decorative aluminium vertical moulded sheet with one side fitted folding foods tablet etc acrylic tray having one box and one drawer fitted with acrylic co lord sheet top with raised edges, drawer and cabinet door front side buffer at rear to avoid collision Base fitted with 4 swivel castors of 100mm non rusted castor wheel with 2 brake,

CARDIAC ICU CHART TROLLEY UMS- Ct03

Overall size top size 900L X 500W stainless steel raised edged fitted top Back side height 1050mm and front side 900mm Machine pressed steel four side raised edged sheet top with second top 12 divided medicine row space third top 8 and fourth top 4 row space Stainless steel Tubular frame mounted on 100mm non rusted castor wheel with 2brake trolley is complete stainless steel Also available pre-treated and epoxy powder coated finish,





DELUX BED SIDE LOCKER WITH ABS TOP UMS-D03

Overall box size : 450L x 450W X 800H mm and drawer size 125HX 300W X 350D mm
 Locker made from machine pressed CRCA ms sheet enclosed on three side which is Pre-treated and epoxy powder coated finis having one box and one drawer fitted with ABS moulded top with raised edges, drawer and cabinet door front side buffer at rear to avoid collision
 Base fitted with 4 swivel castors of 50mm diameter made from synthetic material without brack



DELUX BED SIDE LOCKER WITH SS TOP UMS-P06

Overall box size : 440L x 400W X 800H mm and drawer size 125HX 300W X 350D mm
 Locker made from machine pressed CRCA ms sheet enclosed on three side which is Pre-treated and epoxy powder coated finis having one box and one drawer fitted with SS moulded top with raised edges, drawer and cabinet door front side buffer at rear to avoid collision
 Base fitted with 4 swivel castors of 50mm diameter made Frome synthetic material without brake



STANDARD BED SIDE LOCKER SS TOP UMS-P07

Overall box size : 400L x 400W X 800H mm and cabinet box 300x350 mm
 Locker made from machine pressed CRCA ms sheet locker box with door fitted SS top with raised edges one three sides mild steel tubular lags fitted PVC shoe
 Pre-treated and epoxy powder coated finis,



GEARMATIC HEIGHT ADJUSTBLE OVER BED TABLE UMS-GB03

Overall size 800L X 400W wooden membrane fitted top 800 TO 1100 height adjustment by gear mechanism with handle rectangular tube frame trolley mounted 50 mm synthetic non rusted wheel pre-treated and epoxy powder coated finis,



FIX HEIGHT OVER BED TABLE UMS-GB03

Overall size 800L X 400W MDF board fitted top 1100 fix height Tubular tube frame mounted 50 mm synthetic non rusted wheel pre-treated and epoxy powder coated finis,



DELUX HEIGHT ADJUSTBLE OVER BED TABLE UMS-GA2

Overall size 800L X 400W ABS fitted top 800 TO 1100 height adjustment by pneumatic gas shocker mechanism with handle rectangular tube frame trolley mounted 50 mm synthetic non rusted wheel pre-treated and epoxy powder coated finis,



MANUAL HEIGHT ADJUSTBLE OVER BED TABLE UMS-MO4

Overall size 800L X 400W MDF board fitted top 800 TO 1100 height adjustment by gear mechanism with handle rectangular tube frame trolley mounted 50 mm synthetic non rusted wheel pre-treated and epoxy powder coated finis,