



INSTRUMENT TROLLEY (2) SELVES UMS-517

(A) overall size : 27L X 18W X 36H inch
 (B) overall size: 30L X 24W X 36H inch
 (C) overall size: 42L X 24W X 36H inch
 stainless steel tubular frame mounted on 100 mm non rusted castor wheel with 2 moulded SS self with protective 3 side railings fitted



INSTRUMENT TROLLEY (3) SELVES UMS-518

(A) overall size: 30L X 24W X 36H inch
 (B) overall size: 42L X 24W X 42H inch
 stainless steel tubular frame mounted on 100 mm non rusted castor wheel with 3 moulded SS self with protective 3 side railings fitted



DRESSING TROLLEY (2) SELVES UMS-519

(A) overall size : 27L X 18W X 36H inch
 (B) overall size: 30L X 24W X 36H inch
 stainless steel tubular frame mounted on 100 mm non rusted castor wheel with 2 moulded SS self with protective 3 side railing bowl bucket fitted



OT LINEN TROLLEY UMS-520

(A) overall size : 36L X 18W X 36H inch easy movable handle stainless steel tubular frame mounted on 100 mm non rusted castor wheel with protective 4 side pipe frame fitted canvas cover



MAYO INSTRUMENT TROLLEY DUBLE PILLAR UMS-521

Overall size: tray size 36L X 20W 30 TO 50 inch adjustable height
 SS tubular frame mountad on four non rusted castor wheel with 2 brake



MAYO INSTRUMENT TROLLEY SINGLE PILLAR UM-522

Overall size: tray size 22L X 16W 30 TO 50 inch adjustable height
 SS tubular frame mounted on four non rusted castor wheel with 2 brake