



**QUANTRA DELUX MANUALLY FOWLER BED (2) FUNTION UMS-104**

Technical features		Standard configuration
Bed structure material	MS tube 25x75 pipe 1.4 mm thickness and 1.25 mm thickness CRC sheet with perforated copper wire c02 welding	(1)Head position up down(2) foot position up down
Structure material	MS tube pipe 1.4 mm and 1.25 mm CRC sheet	ABS Moulded panels
Overall bed outer damnation	2000L X980W	Large size2pis tuck away safety side railing
Mattress platform damnation	1900L X 880W	100 mm castor wheel with 2 brake
Height adjustment	Fix height 600 mm	MS CRC perforated sheet
		Hand key with polymer handle
		4 wheel With 2 brake systems
Operating systems	Hand key	4 side bed safety side wafer
Finice epoxy powder coated		4 side iv pole hole
Safe working load	200 kg	Mattress IV stand cylinder trolley optional