



QUANTRA STD MANUALLY (2 FUCTION) FOWLER BED SS PANEL UMS-106

Technical features		Standard configuration
Bed structure material	MS tube 25x75 pipe 1.4 mm thickness and 1.25 mm thickness CRC sheet with perforated copper wire c02 welding	(1)Head position up down (2) foot position up down
Structure material	MS tube pipe 1.4 mm and 1.25 mm CRC sheet	Delux steel panels
Overall bed outer damnation	2000L X980W	Aluminium collapsible safety side railing optional
Mattress platform damnation	1900L X 880W	100 mm castor wheel with 2 brake
Height adjustment	Fix height 600 mm	MS CRC perforated sheet
		Hand key with polymer handle
		4 wheel With 2 brake systems
Operating systems	Hand key	4 side bed safety side wafer
Product finis	Pre-treated and epoxy powder coated	4 side iv pole hole
Safe working load	200 kg	Mattress IV stand cylinder trolley optional