



DELUX BED SIDE LOCKER WITH ABS TOP UMS-D03

Overall box size : 450L x 450W X 800H mm and drawer size 125HX 300W X 350D mm
 Locker made from machine pressed CRCA ms sheet enclosed on three side which is Pre-treated and epoxy powder coated finis having one box and one drawer fitted with ABS moulded top with raised edges, drawer and cabinet door front side buffer at rear to avoid collision
 Base fitted with 4 swivel castors of 50mm diameter made from synthetic material without brack



DELUX BED SIDE LOCKER WITH SS TOP UMS-P06

Overall box size : 440L x 400W X 800H mm and drawer size 125HX 300W X 350D mm
 Locker made from machine pressed CRCA ms sheet enclosed on three side which is Pre-treated and epoxy powder coated finis having one box and one drawer fitted with SS moulded top with raised edges, drawer and cabinet door front side buffer at rear to avoid collision
 Base fitted with 4 swivel castors of 50mm diameter made Frome synthetic material without brake



STANDARD BED SIDE LOCKER SS TOP UMS-P07

Overall box size : 400L x 400W X 800H mm and cabinet box 300x350 mm
 Locker made from machine pressed CRCA ms sheet locker box with door fitted SS top with raised edges one three sides mild steel tubular lags fitted PVC shoe
 Pre-treated and epoxy powder coated finis,