



GEARMATIC HEIGHT ADJUSTBLE OVER BED TABLE UMS-GB03

Overall size 800L X 400W wooden membrane fitted top 800 TO 1100 height adjustment by gear mechanism with handle rectangular tube frame trolley mounted 50 mm synthetic non rusted wheel pre-treated and epoxy powder coated finis,



FIX HEIGHT OVER BED TABLE UMS-GB03

Overall size 800L X 400W MDF board fitted top 1100 fix height Tubular tube frame mounted 50 mm synthetic non rusted wheel pre-treated and epoxy powder coated finis,



DELUX HEIGHT ADJUSTBLE OVER BED TABLE UMS-GA2

Overall size 800L X 400W ABS fitted top 800 TO 1100 height adjustment by pneumatic gas shocker mechanism with handle rectangular tube frame trolley mounted 50 mm synthetic non rusted wheel pre-treated and epoxy powder coated finis,



MANUAL HEIGHT ADJUSTBLE OVER BED TABLE UMS-MO4

Overall size 800L X 400W MDF board fitted top 800 TO 1100 height adjustment by gear mechanism with handle rectangular tube frame trolley mounted 50 mm synthetic non rusted wheel pre-treated and epoxy powder coated finis,