



**Scrub Station Electronic Sensor Operated
(UMS - 509)**



**Scrub Station Foot Operated
(UMS - 510)**

Overall approximate size: 1800L X 600W X900 front side and 1500 back side
Surgical Scrub sinks made of heavy gauge SS 304 material, constructed using ASTM engineering techniques to assure quality and uninterrupted performance. Operation theatre Scrubs are so designed to avoid splashing of water to the neighbouring places and wastage of water. The product ensures trouble free functioning without failure or breakdown. The surgical Scrubs are available in 6 models, in combinations of elbow, knee and Sensor operated, in single, double and three bay options. Water filters, hot and cold mixing options also can be incorporated on demand.

Also available sensor operated

Also available foot operated

Also available hand elbow operated